

# QubeCL System specifications sheet

## 1 Introduction to QubeCL System

### 1.1 Features

- Current noise spectral density  $< 500 \text{ pA}/\sqrt{\text{Hz}}$
- Integrated current noise [10-1000 Hz]  $< 500 \text{ nA}_{\text{RMS}}$
- Current temperature coefficient better than  $5 \text{ ppm}_{\text{FS}}/^\circ\text{C}$
- High current rating up to 2.5 A
- High voltage compliance up to 17.5 V
- Capable of precision modulation from 10 nA
- Fast modulation up to 10 mA @ 2 MHz
- High modulation up to 200 mA @ 600 kHz
- Expansion modules: Temperature stabilization, Lock-in, PLL, Pound-Drever-Hall

### 1.2 Typical applications

- Precision gas spectroscopy
- Locking laser to molecular transition
- Locking laser to high-finesse cavity
- Metrological-grade laser frequency stabilization

### 1.3 General Description

The QubeCL System is a high-performance environment particularly suited to drive state of the art Quantum Cascade laser (QCL) sources, but it can drive any kind of semiconductor laser. The expansion modules allow you to customize the system, obtaining a high-performance instrument in a  $10 \times 10 \text{ cm}^2$  footprint. The system is designed to minimize external wiring and enclose susceptible connections inside its frame, thus reducing noise coupling and maximizing performance.

## 2 Specifications

Interface	Mode	Settings	Notes
USB	Virtual Serial Port	115200 8N1	
WiFi	Access Point/Client	software configurable	

**Table 1:** Communication Interfaces

Parameter	Values	Unit	Notes
Operating temperature	-20 to +40	°C	
Storage temperature	-20 to +70	°C	
Main power supply voltage ( $PS_{LAS}$ )	26	V	No terminal must be connected to ground
TEC power supply voltage ( $PS_{TEC}$ )	26	V	
Monitor power supply voltage ( $PS_{MON}$ )	5.5	V	
Modulation inputs	$\pm 10$	V	

**Table 2:** Absolute Maximum Ratings

### 3 Electrical Specifications

$PS_{LAS} = 24 V_{dc}$ ,  $PS_{TEC} = 12 V_{dc}$ ,  $PS_{MON} = 5 V_{dc}$ , Warm-up time 120 min, room temperature 25 °C, unless otherwise noted.

Parameter	Min	Typ	Max	Unit	Notes
Noise spectral density	100	200	500	pA/ $\sqrt{\text{Hz}}$	test current 200mA
Integrated noise	50	100	200	nA <sub>RMS</sub>	10 Hz - 100 kHz
Integrated noise	250	350	500	nA <sub>RMS</sub>	10 Hz - 1 MHz
Stability (1h)	4	5	6	ppm <sub>FS</sub>	
Stability (12h)	5	6	7	ppm <sub>FS</sub>	
Temperature coefficient	-2	-3	-5	ppm <sub>FS</sub> /°C	
Voltage compliance		17.5		V	

**Table 3:** Current Generator Specs

Parameter	Min	Typ	Max	Unit	Notes
Analog input voltage	0		±10	V	
DC offset	20	30	50	μA	input shorted
Low current					
Noise spect. dens. floor	50		100	pA/ $\sqrt{\text{Hz}}$	input shorted
Output current	0		±5	mA	
Gain		-0.5		mA/V	
Modulation BW (-3 dB)		2		MHz	
High current					
Noise spect. dens. floor	200		2000	pA/ $\sqrt{\text{Hz}}$	input shorted
Output current	0		±100	mA	
Gain		-10		mA/V	
Modulation BW (-3 dB)		600		kHz	

**Table 4:** Analog Current Modulators Specs

Parameter	Min	Typ	Max	Unit	Notes
Waveforms	Sine or Triangle				
Freq. range	$10^{-3}$		$10^6$	Hz	
Min freq. step		250		mHz	
Amplitude control res.		12		bits	
DC offset	20	30	50	$\mu\text{A}$	input shorted
Low current					
Noise spect. dens. floor	50		100	$\text{pA}/\sqrt{\text{Hz}}$	input shorted
Output current	0		$\pm 2.5$	mA	
Modulation BW (-3 dB)		2		MHz	
High current					
Noise spect. dens. floor	200		2000	$\text{pA}/\sqrt{\text{Hz}}$	input shorted
Output current	0		$\pm 75$	mA	
Modulation BW (-3 dB)		600		kHz	

**Table 5:** Mixed Current Generator Specs - Digital

Parameter	Min	Typ	Max	Unit	Notes
Analog input voltage	0		$\pm 2.5$	V	
DC offset	20	30	50	$\mu\text{A}$	input shorted
Low current					
Noise spect. dens. floor	50		100	$\text{pA}/\sqrt{\text{Hz}}$	input shorted
Output current	0		$\pm 2.5$	mA	
Gain		-1		$\text{mA}/\text{V}$	
Modulation BW (-3 dB)		2		MHz	
High current					
Noise spect. dens. floor	200		2000	$\text{pA}/\sqrt{\text{Hz}}$	input shorted
Output current	0		$\pm 75$	mA	
Gain		-30		$\text{mA}/\text{V}$	
Modulation BW (-3 dB)		600		kHz	

**Table 6:** Mixed Current Generator Specs - Analog

Parameter	Min	Typ	Max	Unit	Notes
Thermistor	1		10	k $\Omega$	NTC
TEC current			$\pm 3$	A	
TEC compl. voltage		20		V	
Temp stability		40	400	$\mu$ K	absolute stability in 12h
Temp. coefficient		-90	-100	$\mu$ K/K	

**Table 7:** Temperature Controller Specs

Parameter	Min	Typ	Max	Unit	Notes
RF freq. input	10		250	MHz	
LO freq. input	10		100	MHz	
Input level	-30		0	dBm	
Input stage gain	0		30	dB	
Lock bandwidth			800	kHz	

**Table 8:** PLL Module Specs


Parameter	Min	Typ	Max	Unit	Notes
RF freq. input	1		100	MHz	
LO freq. input	1		100	MHz	
RF Input level	-30		0	dBm	
LO Input level		5		dBm	
Lock bandwidth			800	kHz	

**Table 9:** Pound-Drever-Hall module Specs


Parameter	Min	Typ	Max	Unit	Notes
Modulation frequency		32768		Hz	fixed
Ramp frequency	1		1000	Hz	
Input voltage noise		30		$\mu$ V <sub>RMS</sub>	
Input stage gain	0		60	dB	
Integration time	1		10000	ms	
Lock bandwidth			1	kHz	

**Table 10:** Lock-In module Specs


## 4 Typical performance characteristics




**Figure 1:** Current noise spectral density comparison with low current modulator




**Figure 2:** Current noise spectral density comparison with high current modulator



**Figure 3:** 24h Current stability over ambient temperature



**Figure 4:** 24h Current stability over time



**Figure 5:** Current temperature coefficient



## 5 Current Modulators

The system has two type of current modulators designed to meet the most common needs in terms of modulation. There are two modulation ranges:

- **The low current range** allows to add to the laser even small currents, down to tens of nA, with a high level of precision and very low added noise.
- **The High Current range** allows to add currents up to hundreds of mA to perform large laser frequency scans.

The modulators have a negative gain. If we denote the modulation current with  $I_{mod}$ , the input control voltage  $V_{in}$  and the absolute value of the gain in  $mA/V$  with  $G$ , the modulation current is equal to:

$$I_{mod} = -G \cdot V_{in}.$$

This means that a positive control voltage will subtract current to the laser bias current, and a negative control voltage will add it.

### CAUTION

The high current modulator can provide currents up to 100 mA which can be added or subtracted to the bias current of the laser.

Using negative control voltages the modulation current is added to the bias current, in this case the total current may exceeding the maximum current allowed by the laser. Therefore it is necessary to be very careful to avoid the risk of damaging the laser itself.

Conversely, using positive control voltages it is possible to use the high current modulator safely, because the modulation current is subtracted from the laser bias current.

## 6 Power instructions

### 6.1 Power supply requirements

The system needs three power supplies:

- The power line called  $PS_{LAS}$  has a typical value of 24V and it's used to power the laser driving electronics.
- The power line called  $PS_{TEC}$  has a typical value of 12V and it's used to supply the power circuitry for temperature stabilization.
- The power line called  $PS_{MON}$  has a typical value of 5V and it's used to power the monitor outputs of some boards.

All the voltages used to power the QubeCL **must be floating**, i.e., they **must not be referred to ground**. QubeCL can be powered using common stabilized linear power supplies, the use of switching power supplies can degrade the noise characteristics of the system.

DC input	Min [V]	Typ [V]	Max [V]	Curr. [A]
$PS_{LAS}$	19	24	26	3
$PS_{TEC}$	8	12	26	3
$PS_{MON}$	-	5	5.5	0.5

**Table 11:** Power Supply Reference

#### CAUTION

Be sure to use a FLOATING DC power supply to electrically drive QubeCL system. The power supply terminals must not be connected to ground. Connecting any terminal to ground will damage the instrument!