

# Laser locking with QubeCL and QubeDL

ppqSense Appnote no. 2

**ppqSense s.r.l**

ppqSense Development Team

Revision 1.0 , 2025.05.29



Each QubeCL and QubeDL laser driver may be built in different configurations, allowing it to perform many different operations, from the most basic ones such as current sourcing and temperature stabilization and control, to the most advanced ones.

Particularly, each Qube may be equipped with one locking module, making it capable of actively control and stabilize the wavelength of the driven laser exploiting one of three different locking techniques, one for each available module.

The first available locking module is the **LIA** (Lock-In Amplifier), which enables the Qube to lock a laser to a molecular absorption line.

The **PLL** module allows a Qube to perform an offset-lock to another laser source, used as a reference.

Lastly, the **PDH** module enables a Qube to lock a laser to a high-finesse optical cavity exploiting the Pound-Drever-Hall technique.

This Application Note describes how to use a Qube equipped with the right locking module to lock a laser to the desired reference.

## Contents

<b>1</b>	<b>Frequency lock to a spectroscopic line with LIA module</b>	<b>2</b>
<b>2</b>	<b>Offset phase lock with PLL module</b>	<b>4</b>
<b>3</b>	<b>High finesse cavity lock with PDH module</b>	<b>6</b>
<b>4</b>	<b>Customer support and useful contacts</b>	<b>8</b>

## List of Figures

1	spectroscopic line lock using LIA module . . . . .	2
2	LIA terminals explanation . . . . .	3
3	PLL laser lock setup . . . . .	4
4	PLL terminals explanation . . . . .	5
5	PDH cavity lock setup . . . . .	6
6	PDH terminals explanation . . . . .	7


# 1 Frequency lock to a spectroscopic line with LIA module

An easy way to stabilize the frequency of a laser is to lock its frequency to an atomic or molecular transition. An advantage of this method is that the optical transitions can be quite narrow, typically on the order of a few MHz, however, the thermal Doppler broadening will considerably increase the width of these transitions.

To avoid the Doppler broadening effect, an optical scheme such as the one shown in fig. 1 can be used, saturation spectroscopy allows to achieve a sub-Doppler resolution in order to identify the central frequency of the transition with more precision.

The laser beam to be stabilized is divided into two parts using a beam splitter. Most of the laser radiation is transmitted and is called "pump" beam, a small fraction of the laser radiation is reflected and is called "probe" beam. The two beams have counter-propagating paths and overlap in the cell containing the reference gas. The probe beam is detected after the cell by a photodetector.

This configuration allows to identify the Lamb dip, result of the saturated absorption process, and using a common modulation spectroscopy technique it is possible to obtain the error signal for the laser frequency correction.




**Figure 1:** spectroscopic line lock using LIA module

Using the lock-in module (LIA) installed on the Qube System, it is possible to directly acquire the spectroscopy signal from the photodetector and apply the correction to the laser for frequency locking. In this configuration, the Qube is able to perform internally all the necessary operations to control the locking procedure:

1. Current modulation with frequency and amplitude control
2. Signal demodulation with amplification and integration-time control
3. Processing of the error signal by using a PID controller
4. Addition of the correction current directly to the laser bias current

The Qube Control Software allows to keep under control all the parameters for the optimization of the locking procedure.




**Figure 2:** LIA terminals explanation

## 2 Offset phase lock with PLL module

Frequency locking between two lasers is a technique for transferring the frequency and phase characteristics of a reference laser to a target laser.

It is commonly used in spectroscopic applications and in metrology.

The figure shows a possible optical scheme for an implementation of an offset phase lock between two Quantum Cascade Lasers anyway the Qube System can work with any type of diode laser.




**Figure 3:** PLL laser lock setup

The reference light is superimposed on the target light using a beam splitter.

The overlap of the beams is detected with a wide band photodiode and generates a beat signal at the frequency difference between the target light and the reference, the beat-note frequency is usually in the hundreds of MHz.

The beat-note frequency is a function of the phase difference between the target light and the optical reference and it is constant only if the frequency difference between the target light and the reference light is exactly zero, any phase change will result in a variation of the signal. To perform an offset lock, the beat note is mixed with a local RF reference signal with frequency  $\omega_{offset}$ , at this point the signal obtained from this procedure will be constant only if the frequency difference between the target light and the optical reference is exactly  $\omega_{offset}$ , this is the final error signal used for offset frequency locking.

Using the PLL module integrated in the QuebeCL System, it will be sufficient to send the photodiode signal to the RF IN input and the module will add the correction signal directly to the laser current for frequency/phase stabilization, through the Qube Control Software it is possible to keep under control all the parameters for the optimization of the locking procedure.



**Figure 4:** PLL terminals explanation


### 3 High finesse cavity lock with PDH module

The Pound-Drever-Hall technique is a method to match the emitting optical frequency of a laser to a resonance of a Fabry-Perot cavity. However stable the laser and cavity may be, both the length of the cavity and the frequency of the laser vary over time. PDH locking technique uses the light reflected from the cavity to create an error signal that can be used to correct the laser frequency or change the length of the cavity so that they stay matched and transmission is maximized.

In the PDH technique the laser frequency is modulated, modulated light consists of a carrier frequency and two sidebands. The light reflected from the cavity is measured using a fast photodetector, the reflected signal consists of the two unaltered sidebands together with an out-of-phase carrier component.


The photodetector signal is mixed with a local oscillator, which is in phase with the light modulation. After phase shift and filtering, the resulting electronic signal has a dispersive line shape that can be used as error signal for active stabilization.

For high-finesness cavities, the frequency spacing of the sidebands can be of the order of MHz, in this situation the modulation can be applied directly on the laser current, without the need for external electro-optical modulators.



**Figure 5:** PDH cavity lock setup

The PDH module accepts in input directly the signal coming from the photodiode, processes it and applies the correction to perform the lock directly on the laser current. The Qube Control Software allows to keep under control all the parameters for the optimization of the locking procedure.



**Figure 6:** PDH terminals explanation



## 4 Customer support and useful contacts

The ppqSense staff can be reached for any problem or request concerning the Qube devices. To address any malfunctions or request help while working with our Qube , please write an email to **qube.support@ppqsense.odoo.com**, our staff will reply to the raised ticket as soon as possible.


To obtain quotations, support or updates concerning the purchase of our instruments, please write an email to **info@ppqsense.com**.

All the available documentation concerning both the Qube and all our other products can be found on our website: **www.ppqsense.com**.

The Technical Resources section of the website contains the user manuals, datasheets, brochures and application notes to help users working with the most advanced functionalities of the Qube .

### **ppqSense s.r.l.**

Viale L. Ariosto 492/B  
50019, Sesto Fiorentino (FI)  
Italy  
+39 055 80 23 943



## Revision History

Revision	Date	Author(s)	Description
1.0	2025.05.29	LM	First redaction.